

High-Temperature HPLC Analysis of Cyclosporine on NanoPak-C All-Carbon Columns

Background. Cyclosporine (also known as cyclosporine A or CyA) is an immunosuppressant that has been a drug of choice for the prevention of graft rejection in organ transplant patients since the 1980s [1]. CyA is a cyclic undecapeptide analyzed via RP-HPLC (typically C18 columns) under high-temperature conditions to minimize peak broadening caused by conformational isomerism. Cyclosporine analysis under pharmacopoeial conditions is commonly performed at 80°C, but conventional silica C18 columns can show limited lifetime under these aggressive operating conditions [2], prompting evaluation of alternative stationary phases.

NanoPak-C All-Carbon microbeads are a new class of porous graphitic media that can withstand temperatures up to 200°C [3]. This application note summarizes the evaluation of NanoPak-C All-Carbon columns as an alternative for high-temperature cyclosporine analysis, using European Pharmacopoeia (EP)/ United States Pharmacopoeia (USP)/ Compounded Preparation or Pharmacopoeial (CP)-like solvent systems and conditions.

The goal of the study was to assess whether NanoPak-C All-Carbon could provide suitable chromatographic performance and repeatability at the same elevated temperature while maintaining a retention profile aligned with reference methods.

Table 1. Experimental Conditions	
Cyclosporine solution	1 mg/mL purified cyclosporine in acetonitrile: water (50:50 v/v) solution
Columns	250 mm x 4.6 mm, 6 µm 150 mm x 4.6 mm, 5 µm 50 mm x 4.6 mm, 6 µm
Mobile phase	Phosphoric acid: TBME: Acetonitrile: Water (1: 48: 447: 504, v/v/v/v)
Column	NanoPak-C 250 mm x 4.6 I.D, 6µm
Flow rate	1.9 mL/min
Column oven temp.	80°C
Wavelength	210 nm
Injection volume	20µL
Syringe rinse solvent	Acetonitrile




Figure 1. Representative chromatogram of cyclosporine on (a) 250 X 4.6 mm, (b) 150 X 4.6 mm, (c) 50 X 4.6 mm NanoPak-C All-Carbon columns. (d) Offset chromatograms of seven consecutive injections on the 250 X 4.6 mm column.

Results

On the 250 mm NanoPak-C column, initial retention was 32.980 min at 1.8 mL/min, and after optimization of the acetonitrile/ water composition, the cyclosporine peak eluted at 26.953 min, within the customer's target retention window of 25 to 30 min.

Column stability at 80°C was assessed over seven consecutive injections. Across the sequence,

retention time remained near 27.0 min, tailing factor remained in the range of 1.252 to 1.287, and peak area remained consistent. No significant drift in retention time or peak area was observed across the seven injections, indicating good repeatability and stable performance under the optimized high-temperature conditions.

Taken together, the results reproduced the target pharmacopoeial retention window while maintaining acceptable peak symmetry under 80°C operation.

Shorter columns reduced retention time to 19.081 min on the 150 mm format and 6.387 min on the 50 mm format, supporting their use for faster analyses where monograph-style retention is not required. These formats can therefore be positioned as faster alternatives for in-house or higher-throughput applications. The 250 mm format remains the closest to the reference monograph style.

The data support NanoPak-C All-Carbon HPLC columns as a suitable alternative for routine high-temperature cyclosporine analysis, especially in workflows where reduced lifetime has been observed with conventional silica C18 columns.

Summary

- Matches pharmacopoeial retention window at 80 °C on NanoPak-C All-Carbon 250 × 4.6 mm column.
- Demonstrates stable retention time, peak symmetry, and area over seven injections.
- Shorter column formats enable faster cyclosporine analyses with predictable retention.

Column Selection Guidance

Column format	Observed retention time	Positioning
250 × 4.6 mm, 6 μm	26.95 min	Closest to reference monograph-style method
150 × 4.6 mm, 5 μm	19.08 min	Shorter runtime for faster analysis
50 × 4.6 mm, 6 μm	6.38 min	Very fast screening or in-house high throughput methods

References

- [1] N. Talal, Cyclosporine As An Immunosuppressive Agent For Autoimmune Disease: Theoretical Concepts And Therapeutic Strategies, *Transplant Proc* 20(3 Suppl 4) (1988) 11-5.
- [2] H.A. Claessens, M.A. Van Straten, Review On The Chemical And Thermal Stability Of Stationary Phases For Reversed-Phase Liquid Chromatography, *Journal Of Chromatography A* 1060(1) (2004) 23-41.
- [3] M.J. Parente, B. Sitharaman, Synthesis And Characterization Of Carbon Microbeads, *Acs Omega* 8(37) (2023) 34034-34043.